



MTP Prothesis System



REF CODE	PRODUCT NAME
----------	--------------

METATARSAL COMPONENTS	
-----------------------	--

EXN0001	EXEN MTP METATARSAL COMPONENT SMALL
EXN0002	EXEN MTP METATARSAL COMPONENT MEDIUM
EXN0003	EXEN MTP METATARSAL COMPONENT LARGE

REF CODE	PRODUCT NAME
----------	--------------

PHALANGEAL COMPONENTS	
-----------------------	--

EXN0004	EXEN MTP PHALANGEAL COMPONENT SMALL
EXN0005	EXEN MTP PHALANGEAL COMPONENT MEDIUM
EXN0006	EXEN MTP PHALANGEAL COMPONENT LARGE




INDICATIONS

- Toe arthrosis.

SPECIFICATIONS

- 4 pieces components. Cementless.
- 1 piece metatarsal titanium component, 1 piece phalangeal titanium component, 1 piece titanium metal head, 1 piece polyethylene insert.



REF CODE	PRODUCT NAME
PHALANGEAL COMPONENTS	
EXN0007	EXEN MTP PHALANGEAL INSERT SMALL +0MM
EXN0008	EXEN MTP PHALANGEAL INSERT SMALL +2MM
EXN0009	EXEN MTP PHALANGEAL INSERT SMALL +4MM
EXN0010	EXEN MTP PHALANGEAL INSERT MEDIUM +0MM
EXN0011	EXEN MTP PHALANGEAL INSERT MEDIUM +2MM
EXN0012	EXEN MTP PHALANGEAL INSERT MEDIUM +4MM
EXN0013	EXEN MTP PHALANGEAL INSERT LARGE +0MM
EXN0014	EXEN MTP PHALANGEAL INSERT LARGE +2MM
EXN0015	EXEN MTP PHALANGEAL INSERT LARGE +4MM



REF CODE	PRODUCT NAME
METAL HEAD	
EXN0016	EXEN MTP META HEAD SMALL
EXN0017	EXEN MTP META HEAD MEDIUM
EXN0018	EXEN MTP META HEAD LARGE