


Bio-Force All Inside Meniscal Repair System



REF CODE	PRODUCT NAME
BIO-FORCE MENISCAL REPAIR SYSTEM	
BFP0007	BIO-FORCE PLUS MONO MENISCAL SUTURE STRAIGHT
BFP0008	BIO-FORCE PLUS MONO MENISCAL SUTURE RIGHT
BFP0009	BIO-FORCE PLUS MONO MENISCAL SUTURE LEFT



INDICATIONS

- Meniscal tears repairing.

SPECIFICATIONS

- Attached UHMWPE # 2/0 suture.
- Double PEEK Implants.
- Adjustable self locking knotted suture system.