

Peek-Force Plus Screw Type Suture Anchors



REF CODE	PRODUCT NAME
PEEK-FORCE PLUS SCREW TYPE SUTURE ANCHORS	
PFP0011	PEEK-FORCE PLUS SCREW TYPE SUTURE ANCHOR 5.5 MM
PFP0013	PEEK-FORCE PLUS SCREW TYPE SUTURE ANCHOR WITH NEEDLE 5.5 MM
PFP0010	PEEK-FORCE PLUS SCREW TYPE SUTURE ANCHOR 3.5 MM
PFP0012	PEEK-FORCE PLUS SCREW TYPE SUTURE ANCHOR WITH NEEDLE 3.5 MM



INDICATIONS

- 3.5mm anchors for PASTA lesion repair, Bankart repair and SLAP surgery
- 5.5mm anchors for Rotator Cuff repair surgery

SPECIFICATIONS

- 3.5mm, 5.5mm anchors are loaded; 2 x UHMWPE sutures 900mm length USP # 2