



REF CODE	PRODUCT NAME
TI-FORCE PLUS LIFT DOUBLE BUTTON	
TFP0023	Ti-FORCE PLUS LIFT DOUBLE BUTTON



INDICATIONS

- Ankle syndesmosis repairs

SPECIFICATIONS

- Oval titanium button with four-holes, round titanium button with two-holes.
- The main body of the implants consist of #5 USP UHMWPE sutures.
- There is also USP #2 UHMWPE guidance suture.
- The round button of the implants are fully compatible to the screw holes of all brands of trauma plates.
- Besides, It is suitable for use with or without trauma plates according to the surgical application techniques.
- They have a self-locking loop suture design.