

## Ti-Force ACL Fixation Devices



| REF CODE                             | PRODUCT NAME                       |
|--------------------------------------|------------------------------------|
| <b>Ti-FORCE ACL FIXATION DEVICES</b> |                                    |
| TFR0026                              | Ti-FORCE ACL FIXATION DEVICE 15 MM |
| TFR0027                              | Ti-FORCE ACL FIXATION DEVICE 20 MM |
| TFR0028                              | Ti-FORCE ACL FIXATION DEVICE 25 MM |
| TFR0029                              | Ti-FORCE ACL FIXATION DEVICE 30 MM |
| TFR0030                              | Ti-FORCE ACL FIXATION DEVICE 35 MM |
| TFR0031                              | Ti-FORCE ACL FIXATION DEVICE 40 MM |



### INDICATIONS

- ACL-PCL surgery.

### SPECIFICATIONS

- The Implant has single titanium button, pre-sized loaded UHMWPE continious loop.
- For guidance 2 x UHMWPE sutures 900mm lenght USP #2