

Crosspes Screws



REF CODE	PRODUCT NAME
CROSSPES SCREWS	
CRSP70	CROSSPES SCREW Ø7MM
CRSP80	CROSSPES SCREW Ø8MM
CRSP90	CROSSPES SCREW Ø9MM
CRSP100	CROSSPES SCREW Ø10MM



INDICATIONS

- Severe pronation
- Calcaneal valgus deformity
- Talus in plantarflexion
- Failed correction with long-term orthotic treatment
- Congenital and painful flat foot deformity
- Treatment of tarsal unions
- Subtalar instability
- Posterior tibial tendon dysfunction
- Paralytic flat foot deformity

SPECIFICATIONS

- Titanium, cannulated.
- Inside holes for osteointegration.
- Talus designed to prevent normal subtalar joint movement.