



REF CODE	PRODUCT NAME
LUMBAR PEEK TLIF CAGES	
AVS2426	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H7
AVS2427	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H8
AVS2428	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H9
AVS2429	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H10
AVS2430	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H11
AVS2431	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H12
AVS2432	AVES PLUS LUMBAR PEEK CAGE TLIF 0° H13

INDICATIONS

- In disc herniation
- In unstable discopathy degeneration
- In restoring disk height

SPECIFICATION

- Bone retention / support area is maximized before fusion thanks to titanium pins and special surface design
- PEEK material
- Graft spaces are optimal