

Force Halluxvalgus Wedge Plate



REF CODE	PRODUCT NAME
FORCE HALLUX VALGUS WEDGE PLATE	
FRC2006	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE RIGHT 0
FRC2007	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE RIGHT 2
FRC2008	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE RIGHT 4
FRC2009	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE RIGHT 6
FRC2010	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE LEFT 0
FRC2011	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE LEFT 2
FRC2012	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE LEFT 4
FRC2013	FORCE HALLUX VALGUS WEDGE OSTEOTOMY PLATE LEFT 6



INDICATIONS

- Hallux Valgus deformities.

SPECIFICATIONS

- 1.2 mm thick plates and self tapping screws.
- The plates are bended according to the bone's natural shapes.
- Locked design.