

Force Hallux Valgus Chevron Osteotomy Plate



REF CODE	PRODUCT NAME
FORCE HALLUX VALGUS ANGLED CHEVRON PLATE	
FRC2000	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE RIGHT 33
FRC2001	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE RIGHT 35
FRC2002	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE RIGHT 37
FRC2003	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE LEFT 33
FRC2004	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE LEFT 35
FRC2005	FORCE HALLUX VALGUS CHEVRON ANGLED OSTEOTOMY PLATE LEFT 37



INDICATIONS

- Hallux Valgus deformities.

SPECIFICATIONS

- 1.5 mm thick plates and self tapping screws.
- The plates are bended according to the bone's natural shapes.
- Locked design.