

Knotless Type

Peek-Force Knotless Peek-Ti Anchors



REF CODE	PRODUCT NAME
PEEK-FORCE KNOTLESS PEEK-Ti ANCHORS	
PFR0034	PEEK-FORCE KNOTLESS PEEK-Ti ANCHOR 2.8 MM
PFR0035	PEEK-FORCE KNOTLESS PEEK-Ti ANCHOR 3.5 MM
PFR0037	PEEK-FORCE KNOTLESS PEEK-Ti ANCHOR 5.0 MM



INDICATIONS

- 2.8mm, 3.5mm anchors for Bankart and SLAP surgery
- 5.0mm anchors for Rotator Cuff repair surgery

SPECIFICATIONS

- Anchor bodies are made of peek material. The tip of the anchors are made of Titanium.
- Anchors do not contain sutures. It is recommended to be used externally with ALEDA brand USP #2 UHMWPE sutures.
- Sliding design.