

Knee Arthroscopic Implant Types and Ancillary Products

Force Plus Lift up ACL-PCL Fixation Device with Loop



REF CODE	PRODUCT NAME
	FORCE PLUS LIFT UP ACL-PCL FIXATION DEVICE WITH LOOP
TFR0036	FORCE PLUS LIFT UP ACL-PCL FIXATION DEVICE WITH LOOP



INDICATIONS

- ACL-PCL surgery.

SPECIFICATIONS

- The Implant has single titanium button and loaded; for loop 1 x UHMWPE sutures USP #7.
- For guidance 2 x UHMWPE sutures 900mm length USP #2.
- Adjustable from upper side and self locking system.