

## Force Lift down ACL-PCL Fixation Device with Loop



REF CODE	PRODUCT NAME
TFR0037	FORCE LIFT DOWN ACL-PCL FIXATION DEVICE WITH LOOP



### INDICATIONS

- ACL-PCL surgery

### SPECIFICATIONS

- The Implant has single titanium button and loaded; for loop 1 x UHMWPE sutures USP #7.
- For guidance 2 x UHMWPE sutures 900mm length USP #2.
- Adjustable from lower side and self locking system.