

# FORCE LIFTDOWN ACL/PCL FIXATION DEVICE WITH LOOP



This products are manufacturing by Titanium alloy and UHMWPE fibers for ACL and PCL repair. Button measurements are; lenght 12mm widht 3,90mm and thickness 1.5mm. They can be use with all brand acl instruments with 4.5mm or 5mm cannulated drills. Loop lenghts are adjustable that what the surgeons needs. Loop ropes width is USP #7, and adjusting ropes are USP # 2. This implant's loop lenghts adjusting and locking from tibial tunnel or arthroscopic portals of the knee by way of special self locking design with free rope ends.

*"POWERFUL, FAST AND RELIABLE"*

[www.aledamakina.com](http://www.aledamakina.com)



**REF CODE**

**PRODUCT NAME**

**FORCE LIFTDOWN ACL/PCL FIXATION DEVICE WITH LOOP**

YN00235

FORCE LIFTDOWN ACL/PCL FIXATION DEVICE WITH LOOP

Çamlıca Mah. Anadolu Bulvarı R Blok Apt. No: 20 R/6  
Yenimahalle | ANKARA | TURKEY

**Phone:** +90 312 397 70 41 | **Fax:** +90 312 397 70 42

info@aledamakina.com