

Knee Arthroscopic Implant Types and Ancillary Products

Fireforce Monopolar Radiofrequency (RF) Ablation Probes



REF CODE	PRODUCT NAME
FIREFORCE MONOPOLAR RADIOFREQUENCY (RF) ABLATION PROBES	
MRF45AP	FIREFORCE MONOPOLAR RADIOFREQUENCY (RF) ABLATION PROBES 45°
MRF90AP	FIREFORCE MONOPOLAR RADIOFREQUENCY (RF) ABLATION PROBES 90°
MRF180AP	FIREFORCE MONOPOLAR RADIOFREQUENCY (RF) ABLATION PROBES 180°

INDICATIONS

- For Arthroscopic joint surgery; (Knee joint, Shoulder joint, Ankle joint) for cutting, shrinking, coagulations.

SPECIFICATIONS

- Monopolar with ceramic ring, probe's body coated stainless steel.
- Ø3.0mm x 165mm
- These products; It works compatible with all brands of cautery generators.

