

# Knotless Type <

## > Peek-Force Knotless Peek-Ti Anchors



REF CODE	PRODUCT NAME
<b>PEEK-FORCE KNOTLESS PEEK-Ti ANCHORS</b>	
PFR0036	PEEK-FORCE KNOTLESS PEEK-Ti ANCHOR 4.5 MM
PFR0038	PEEK-FORCE KNOTLESS PEEK-Ti ANCHOR 5.5 MM



### INDICATIONS

- 4.5mm, 5.5mm anchors for Rotator Cuff repair surgery

### SPECIFICATIONS

- Anchor bodies are made of PEEK. The tip of Anchor bodies are made of Titanium.
- Anchors do not contain sutures. It is recommended to be used externally with ALEDA brand USP #2 UHMWPE sutures.
- No need instruments for drilling.