

PEEK-FORCE KNOTLESS PEKK ANCHORS



REF CODE	PRODUCT NAME
PEEK-FORCE SUTURE ANCHORS	
PFR0003	PEEK-FORCE KNOTLESS PEKK ANCHOR Ø3.5 MM
PFR0006	PEEK-FORCE KNOTLESS PEKK ANCHOR Ø5.5 MM
PFR0012	PEEK-FORCE KNOTLESS PEKK ANCHOR WITH NEEDLE Ø3.5MM
PFR0013	PEEK-FORCE KNOTLESS PEKK ANCHOR WITH NEEDLE Ø5.5MM



INDICATIONS

- 3.5mm anchors for PASTA lesion repair, Bankart repair and SLAP surgery
- 5.5mm anchors for Rotator Cuff repair surgery

SPECIFICATIONS

- 3.5mm, 5.5mm anchors are loaded; 2 x UHMWPE sutures 900mm length USP # 2