

Knotless Type

Anchors

Peek-Force Plus Knotless Peek Anchors



REF CODE	PRODUCT NAME
PEEK-FORCE PLUS KNOTLESS PEEK ANCHORS	
PFP0002	PEEK-FORCE PLUS KNOTLESS PEEK ANCHOR 2.8 MM
PFP0003	PEEK-FORCE PLUS KNOTLESS PEEK ANCHOR 3.5 MM
PFP0004	PEEK-FORCE PLUS KNOTLESS PEEK ANCHOR 4.5 MM
PFP0005	PEEK-FORCE PLUS KNOTLESS PEEK ANCHOR 5.0 MM
PFP0006	PEEK-FORCE PLUS KNOTLESS PEEK ANCHOR 5.5 MM



INDICATIONS

- 2.8mm anchors for PASTA lesion surgery, Bankart and SLAP surgery
- 3.5mm anchors for Bankart and SLAP surgery
- 4.5mm, 5.0mm, 5.5mm anchors for Rotator Cuff repair surgery

SPECIFICATIONS

- Anchor bodies are completely composed of peek material.
- Anchors do not contain sutures. It is recommended to be used externally with ALEDA brand USP #2 UHMWPE sutures.