


Knotless Type <

Anchors

> Talon Peek-Ti Knotless Anchor



REF CODE	PRODUCT NAME
TALON PEEK-Ti KNOTLESS ANCHOR	
TLN0001	TALON PEEK-Ti KNOTLESS ANCHOR
SIZE	5.5 MM



INDICATIONS

- 5.5mm Talon anchors for Rotator Cuff repair surgery
- It is for especially osteoporotic patients.

SPECIFICATIONS

- 5.5mm Talon anchors; 2 x UHMWPE sutures 900mm length USP # 2
- Anchors do not contain sutures. It is recommended to be used externally with ALEDA brand USP #2 UHMWPE sutures.