

## ➤ Ti-Force Full Thread Anchors



REF CODE	PRODUCT NAME
<b>Ti-FORCE FULL THREAD ANCHORS</b>	
TFR0043	Ti-FORCE FULL THREAD ANCHOR 4.5 MM
TFR0044	Ti-FORCE FULL THREAD ANCHOR WITH NEEDLE 4.5 MM
TFR0045	Ti-FORCE FULL THREAD ANCHOR 5.5 MM
TFR0046	Ti-FORCE FULL THREAD ANCHOR WITH NEEDLE 5.5 MM



### INDICATIONS

- 4.5mm - 5.5mm anchors for Rotator Cuff repair surgery

### SPECIFICATIONS

- 4.5mm, 5.5mm anchors are loaded; 2 x UHMWPE sutures 900mm length USP # 2