

## Exen Titanium Locking Cerclage Elastic Cable System



### INDICATIONS

- For Upper and lower extremities fractures and general orthopedic trauma surgery, especially Femur (including periprosthetic fractures), prophylactic banding during total hip joint procedures.

### SPECIFICATIONS

- Titanium body with excellent safe locking system due to unique design.
- Dimensions; height: 6mm, width: 11mm, length: 13mm.
- 1.50mm x 1000mm with core braided UHMWPE sutures.
- Suitable with all brands trauma plates and their transverse cannals.

REF CODE	PRODUCT NAME
<b>EXEN TITANIUM LOCKING CERCLAGE ELASTIC CABLE SYSTEM</b>	
EXN1714	EXEN TITANIUM LOCKING CERCLAGE ELASTIC CABLE SYSTEM